Epic™

OVER-THE-WIRE

Self-Expanding Nitinol Vascular Stent with Delivery System

Directions for Use	2
Instrucciones de uso	7
Mode d'emploi	13
Gebrauchsanweisung	19
Istruzioni per l'uso	25
Gebruiksaanwijzing	31
取扱説明書	37

TABLE OF CONTENTS

WARNING	3
DEVICE DESCRIPTION	3
INTENDED USE	3
CONTRAINDICATIONS	3
WARNINGS	3
PRECAUTIONS	3
MAGNETIC RESONANCE IMAGING (MRI)	4
ADVERSE EVENTS	4
HOW SUPPLIED	5
RECOMMENDED MATERIALS	5
OPERATIONAL INSTRUCTIONS	5
Patient Preparation	5
Inject Contrast Media	5
Evaluate and Mark the Stenosis	5
Select Proper Stent System	5
Preparation of Stent Delivery System	5
Figure 1	5
Delivery Procedures.....	5
Stent Deployment Procedure (Reference Figure 1).....	6
Post Procedure	6
REFERENCES	6
WARRANTY	6

Epic™

OVER-THE-WIRE

Self-Expanding Nitinol Vascular Stent with Delivery System

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

The Epic Nitinol Vascular Stent System is comprised of two components: the implantable endoprosthesis and the stent delivery system. The stent is a laser cut self-expanding stent composed of a nickel titanium alloy (nitinol). On both the proximal and distal ends of the stent, radiopaque markers made of tantalum increase visibility of the stent to aid in placement. The stent is constrained within a 6F 2.1 mm maximum OD delivery system. The delivery system is a coaxial design with an exterior shaft to protect and constrain the stent prior to deployment. The delivery system is compatible with 0.035 in (0.89 mm) guidewires.

When ready to be implanted, the stent is deployed by retracting the exterior shaft of the delivery system. A radiopaque marker at the distal end of the delivery system aids in visibility during deployment. As the stent is exposed to body temperature, it expands to appose the vessel wall.

The Epic Nitinol Vascular Stent System is available in a variety of stent diameters and lengths. The delivery system is also offered in two shaft lengths.

Please see the product label for the specific delivery system length, stent diameter and stent length.

This product contains no detectable latex.

INTENDED USE

The Epic Nitinol Vascular Stent System is indicated for the treatment of peripheral vascular lesions and obstructions.

CONTRAINDICATIONS

Generally, contraindications for percutaneous transluminal angioplasty (PTA) are also contraindications for stent placement. Contraindications for the use of the Epic Nitinol Vascular Stent System include, but may not be limited to:

- Patients with highly calcified lesions resistant to PTA.
- Patients with a target lesion with a large amount of adjacent acute or subacute thrombus.
- Patients with uncorrected bleeding disorders or patients who cannot receive anticoagulation or antiplatelet aggregation therapy.
- Patients with perforated vessels evidenced by extravasation of contrast media.
- A lesion that is within or adjacent to an aneurysm.
- Patients with excessive vessel tortuosity.

WARNINGS

Stenting across a bifurcation or side branch could compromise future diagnostic or therapeutic procedures.

Persons allergic to nickel may suffer an allergic response to this implant.

The stent is not designed for repositioning.

Once the stent is partially deployed, it cannot be “recaptured” or “resheathed” using the stent delivery system.

As with any type of intravascular implant, infection secondary to contamination of the stent, may lead to thrombosis, pseudoaneurysm, or rupture into a neighboring organ or into the retroperitoneum.

The stent may cause thrombus or distal emboli to migrate from the site of the implant down the arterial lumen.

Do not use if the exposure indicator dot is red, indicating that the stent expansion may have been compromised.

Do not use if the temperature exposure indicator dot is missing.

Do not expose to organic solvents (e.g. alcohol).

PRECAUTIONS

- Use prior to the ‘use by’ date on the package label.
- The device is intended for use by physicians who have received appropriate training.
- The delivery system is not designed for use with power injection systems.
- Do not use a kinked delivery system.
- Only advance the stent delivery system over a guidewire.
- When catheters are in the body, they should be manipulated only under fluoroscopy. Radiographic equipment that provides high quality images is needed.
- When treating multiple lesions, the most distal lesions should be stented first followed by the stenting of proximal lesions. Stenting in this order eliminates the need to cross and reduces the chance of dislodging stents that have already been placed.

- When multiple stents are required, if placement results in metal to metal contact, stent materials should be of similar composition and overlapped at least 5 mm.
- The stent delivery system is not intended for arterial blood monitoring.
- Prior to completion of the procedure, utilize fluoroscopy to ensure proper positioning of the stent. If the target lesion is not fully covered, use additional stents as necessary to adequately treat the lesion.
- The minimally acceptable sheath French size is printed on the package label. Do not attempt to pass the stent delivery system through a smaller size sheath introducer than indicated on the label.
- In the event of thrombosis of the expanded stent, thrombolysis and PTA should be attempted.
- In the event of complications such as infection, pseudoaneurysms or fistulization, surgical removal of the stent may be required. Standard surgical procedure is appropriate.
- Recrossing a stent with adjunct devices must be performed with caution.
- Premature removal of the safety lock may result in an unintended deployment of the stent.

MAGNETIC RESONANCE IMAGING (MRI)

Through non-clinical testing, the Epic™ Stent has been shown to be MR Conditional, at field strengths from 1.0 to 3.0 Tesla, with maximum whole body averaged specific absorption rates (SAR) and maximum continuous MRI scan times stated in the table below. Other conditions of the testing include the following:

- Spatial gradient fields of 70 mT/cm or less
- A rate of change of magnetic field (dB/dt) of 60 T/s or less

The following table lists the conditions that produce a maximum temperature rise of 2.0 degrees C.

Maximum continuous MRI scan time based on 1.0 to 1.4 W/Kg SAR Level

Total Stent Length (including overlap) (mm)	Static Magnetic Field (Tesla)	SAR (W/Kg)	Maximum Continuous MRI Scan Time
120	1.0	1.4	15 min.
	1.5	1.0	10 min. 45 sec.
	3.0	1.1	15 min.
180	1.0	1.0	1 min. 45 sec.
	1.5	1.0	3 min.
230	1.0	1.0	1 min. 15 sec.
	1.5	1.0	2 min. 30 sec.
	3.0	1.1	15 min.

The RF heating testing was conducted on a Philips NT Intera™ 1.0 Tesla scanner, a Philips Intera 1.5 Tesla scanner, and a Siemens Magnetom Trio™ 3 Tesla scanner. RF heating with scanners other than those listed is unknown.

The effect of heating in the MRI environment for stents with simulated fractures has been tested and found to be similar to single stents.

MRI at 3T or less may be performed immediately following the implantation of the Epic Stent. The Epic Stent should not migrate in this MRI environment. This stent has not been evaluated to determine if it is MR Conditional beyond these conditions.

MR imaging quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the stent.

ADVERSE EVENTS

Procedures requiring percutaneous catheter introduction should not be attempted by physicians unfamiliar with the possible complications. Complications may occur at any time during or after the procedure. Potential complications may include, but are not limited to:

- Abscess
- Acute/subacute stent thrombosis
- Allergic reaction (to drug, contrast, device, or other)
- Amputation
- Aneurysm
- Angina/coronary ischemia
- Arrhythmia
- AV fistula
- Death
- Drug reactions
- Embolization (air, plaque, thrombus, device, or other)
- Entanglement of delivery system in deployed stent
- GI bleeding from anticoagulant/antiplatelet medication
- Hemorrhage/hematoma
- Hypotension/hypertension
- Myocardial Infarction (MI)
- Need for urgent intervention or surgery
- Pseudoaneurysm formation
- Renal insufficiency or failure
- Restenosis of stented artery
- Rupture of the retroperitoneum or neighboring organ
- Sepsis/infection
- Stent fracture
- Stent migration
- Stent misplacement/jumping
- Stroke
- Thrombus
- Tissue necrosis/ischemia
- Vasospasm
- Vessel injury (perforation, dissection, intimal tear, rupture, or other)
- Vessel occlusion

HOW SUPPLIED

Handling and Storage

Store in a cool, dry dark place and do not exceed 55°C.

Do not use if package is opened or damaged.

RECOMMENDED MATERIALS

- Micropuncture kit
- 0.035 in (0.89 mm) guidewire of appropriate length
- Introducer sheath of appropriate size and length and equipped with a hemostatic valve
- Syringe (10 ml (cc) for prepping the stent delivery system)

OPERATIONAL INSTRUCTIONS

Patient Preparation

The percutaneous placement of a vascular self-expanding stent in a peripheral vascular stenotic or obstructed artery should be done in an angiography procedure room equipped with the appropriate imaging equipment. Patient preparation and sterile precautions should be the same as for any angioplasty procedure. Angiography should be performed to map out the extent of the lesion(s) and the collateral flow. Access vessels must be sufficiently patent, or sufficiently recanalized, to proceed with further intervention. If thrombus is present or suspected, thrombolysis should precede stent deployment using standard acceptable practice.

Inject Contrast Media

Perform angiogram using standard technique.

Evaluate and Mark the Stenosis

Observe fluoroscopically the most distal view of the stenotic or obstructed artery.

Select Proper Stent System

1. Measure the diameter of the reference vessel (proximal and distal to the lesion or obstruction). Select a stent per the table below to achieve a secure placement:

Reference Vessel Diameter (mm)	Self Expanding Stent Diameter (mm)
4.0	5.0
4.0-5.0	6.0
5.0-6.0	7.0
6.0-7.0	8.0
7.0-8.0	9.0
8.0-9.0	10.0
9.0-11.0	12.0
11.0-13.0	14.0

2. Measure the entire length of the actual lesion and select the proper length of the stent(s) to be deployed. To help ensure adequate apposition, it is recommended that the length of the stent be chosen so that the ends of the stent extend at least 5mm beyond both ends of the lesion into healthy tissue. Should more than one stent be required to cover the lesion, place the distal stent first followed by the proximal stent(s), allowing for at least 5 mm of stent overlap.

3. Estimate the distance between the lesion and the entry site to select the proper stent delivery system length.

Preparation of Stent Delivery System

1. Open the outer box to reveal the pouch containing the stent delivery system.
2. Check the temperature exposure indicator on the pouch label to confirm that the product has not been compromised. See Warnings section.
3. After careful inspection of the pouch looking for damage to the sterile barrier, carefully peel open the pouch and extract the stent delivery system tray.
4. Open the door on the tray that contains the handle.
5. Carefully withdraw the stent delivery system from the tray by grasping the handle of the delivery system.
6. Examine the device for any damage. If it is suspected that the sterility or performance of the device has been compromised, the device should not be used.
7. If the safety lock ③ (Reference Figure 1) is not attached to the device, verify that the stent is fully constrained in the delivery system and place safety lock in position as shown in Figure 1.
8. Attach a 10 ml (cc) syringe filled with saline to the luer ⑦ (Reference Figure 1) on the handle. Apply positive pressure. Continue to flush until saline appears at the distal end of both the guidewire lumen and the sheath–tip junction. Note: If the luer is not present when received, the device may require greater than 10cc to flush entire device.
9. Remove the flushing luer ⑦.


Figure 1.

Delivery Procedures

1. Gain arterial access utilizing a 6F (2.1mm) sheath with a hemostatic valve.

Precaution: Always use an introducer or guide sheath for the implant procedure, to protect the access site. An introducer sheath of a 6F or larger size is recommended.

2. Pass a 0.035 in (0.89 mm) guidewire of appropriate length across the target lesion or obstruction.
3. Pre-dilate the lesion as necessary with a balloon dilatation catheter using conventional technique.

Precaution: Physicians should use judgement based on experience in dilating arterial lesions and/or obstructions. Never force a balloon catheter to inflate to the point of risking dissection of the arterial wall.

After the lesion has been properly dilated, remove the dilatation catheter, leaving the guidewire with the tip distal to the lesion for stent system advancement.

4. Place the Epic™ Nitinol Vascular Stent System over the guidewire. Advance the delivery system as a unit through the hemostatic valve of the introducer sheath.

Stent Deployment Procedure (Reference Figure 1)

1. Advance the delivery system until the stent radiopaque markers ① and ② are centered over the target lesion.

Precaution: If unable to initiate release of the stent or if strong resistance is met with the introduction of the delivery system, remove the entire system from the patient and introduce a new system. Do not use a power injector through the delivery system for angiography.

2. Remove the safety lock positioned on the rack ③ by pulling vertically. Confirm that the radiopaque markers are still properly positioned across the target lesion. Keep the entire length of the delivery system as straight as possible and maintain slight backward tension on the delivery system during deployment. Note: If repositioning of the delivery system is required, reinserting the safety lock will prevent inadvertent deployment.

Deploy the stent by one of the following methods:

- A. Roll the thumb wheel ⑤ of the deployment handle in a proximal direction. Continue to roll thumb wheel until the radiopaque marker of the exterior shaft ⑥ passes the proximal radiopaque markers of the stent resulting in full deployment. Do not continue to roll the thumb wheel after the stent is fully deployed.
 - B. Grasp the manual pull grip ④ and pull toward the deployment handle. Continue to pull back until the radiopaque marker of the exterior shaft ⑥ passes the proximal radiopaque markers of the stent resulting in full deployment. Do not continue to pull back the manual grip after the stent is fully deployed. Note: Do not restrict movement of the thumb wheel ⑤ otherwise deployment difficulties could be encountered.
 - C. (Combination) Any combination of A and B can be used to achieve full deployment.
3. When released from the delivery system, the stent will immediately expand into position. View the delivery system under fluoroscopy, ensuring that the exterior shaft marker band ⑥ has crossed the proximal stent markers. The delivery system can now be withdrawn. Use caution when withdrawing the stent delivery system and always manipulate under fluoroscopy. If unusual resistance is felt, readvance and rotate the delivery system in an attempt to center the delivery system within the vessel, then carefully attempt repeat withdrawal.

4. If incomplete expansion exists within the stent at any point along the lesion, balloon dilatation can be performed utilizing standard PTA technique.

Precaution: Never dilate the stent using a balloon that is larger in diameter than the unconstrained diameter of the stent.

5. Withdraw guidewire and sheath from patient and establish hemostasis per conventional technique.

Post Procedure

Assess patient for hematoma and/or other signs of bleeding at the puncture site.

REFERENCES

The physician should consult recent literature on current medical practice on stent implantation.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

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© 2008 Boston Scientific Corporation or its affiliates. All rights reserved.
Boston Scientific International S.A.